

For the use only of a Registered Veterinary Practitioner or a Hospital or a Laboratory or a Farm.

Lactic Acid and Salicylic Acid Ear Cleanser

## CL-Ear™ Cleanser

Ear Cleanser for  
Dog and Cats

### COMPOSITION:

Each ml contains:  
Lactic Acid USP 2.5% w/v  
Salicylic Acid USP 0.1% w/v  
Excipients q.s.

### DESCRIPTION:

**CL-Ear Cleanser** is an ear cleanser containing kerolytic agent salicylic acid and lactic acid.

For auricular use in Dogs & Cats only.

### INDICATIONS:

**CL-Ear Cleanser** is formulated to remove crust and excessive cerumen (Ear Wax). For routine ear cleansing in healthy ear canals or prior to administration of other ear preparations in the treatment of otitis externa.

### DOSAGE AND ADMINISTRATION:

#### Dosage:

From two to three times a week to daily for routine use, before each ear treatment in the case of otitis externa, or as directed by your veterinarian.

#### Administration:

For auricular use/ To be applied directly into the ear canals.

#### Directions for use:

Remove the existing cap and seal and use twist-on cap for application. Turn the cap nozzle to open the product and then apply liberally into the ear canal. Gently rub the base of the ear and then wipe the interior of the ear flap with cotton or cloth moistened with ear cleanser.

### PHARMACODYNAMIC PROPERTIES:

**Salicylic acid and Lactic acid:** They are kerolytic agent with antibacterial properties. Salicylic acid is believed to have the property of solubilising the cell surface proteins, resulting in desquamation of keratolytic debris.

**Lactic acid:** It may diminish corneocyte adhesion by interfering with the formation of ionic bonds.

### WARNING & PRECAUTIONS:

Special precautions to be taken by the person administering the veterinary medicinal product to animals. Avoid contact with the product. Do not touch the dropper tip to any surface since it may contaminate the solution. Wash hands carefully after applying the product. In case of accidental contact with the eyes, rinse with copious amounts of water. Do not use in ears with ruptured tympanic membranes (eardrums).

### CONTRAINDICATIONS:

Do not administer to animals other than those for which it is indicated.  
Hypersensitivity to any of the Ingredients.

### STORAGE:

Store in a cool and dry place.

### PRESENTATION:

- 1) 50 ml HDPE Bottle packed in a Primary Carton along with the Pack Insert and a twist dropper cap.
- 2) 100 ml HDPE Bottle packed in a Primary Carton along with the Pack Insert and a twist dropper cap.

**VETERINARY. NOT FOR HUMAN USE.  
FOR ANIMAL TREATMENT ONLY.  
FOR EXTERNAL APPLICATION ONLY.  
NOT FOR INJECTION.**

**SHAKE WELL BEFORE USE.**

**KEEP OUT OF REACH OF CHILDREN & PETS, AWAY FROM FOOD.**

Carefully read the accompanying instructions before use.

# Veko®

Manufactured & Marketed by:

**Veko Care Pvt. Ltd.**

Plot No. E-48 & 49, MIDC, Ranjangaon, Dist-Pune.

Maharashtra, India. Pin Code: 412220

Customer care: care@vekocare.com

TM : Trademark of Veko care



VPP1E002